

# CARE PUMP




## **LITE8** *Economic line*

<b>Number of chambers</b>	8
<b>Max. pressure [mmHg]</b>	220
<b>Number of proprietary therapeutic pre-programs of the medical organization</b>	8
<b>Precise regulation of pressure and time</b>	in steps of 1
<b>Number of independent channels</b>	1
<b>Adjustment of the speed of the treatment</b>	1-5
<b>Cuffs</b>	system overlapping
<b>Display</b>	5" TN Icd color touch screen
<b>Weight [kg]</b>	3,2
<b>Dimensions [mm]</b>	300x260x160
<b>Power supply</b>	100÷240 V, 50÷60 Hz
<b>Power consumption</b>	40-70 W
<b>Warranty</b>	1 year

# Accessories




## Included:




power supply

## Optional:




8-chamber arm cuff

CPARM8T




8-chamber waist cuff

CPWAIST8T



8-chamber left leg cuff

CPLEG8TL




8-chamber right leg cuff

CPLEG8TR




8-chamber short pants

CPSHORTPANTS8T




Extension zipper

CPLEGEXT8




Double splitter

CPDOUBLECONNECTOR8




Patient's safety button

SAFETY BUTTON




8-chamber double cuff for arms, shoulders and chest

CPDOUBLEARM8T




8-chamber left arm, shoulder and chest cuff

CPARMSB8TL




8-chamber right arm, shoulder and chest cuff

CPARMSB8TR




8-chamber full pants

CPLEG8TL




2 x Extension zippers for 8-chamber double cuff for arms, shoulders and chest

CPDOUBLEARMEXT8




2 x Extension zippers for 8-chamber short pants

CPSHORTPANTSEXT8



2 x Extension zippers for 8-chamber full pants

CPLONGPANTSEXT8



Double splitter (long - 150 cm) for 8-chamber devices

CPDOUBLECONNECTOR8L

# Therapeutic pre-programs

of the medical organization



mode **A**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	50 min.
<b>Hold</b>	3 s
<b>Interval</b>	3 s

mode **B**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	50 min.
<b>Hold</b>	3 s
<b>Interval</b>	3 s

mode **C**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	50 min.
<b>Hold</b>	3 s
<b>Interval</b>	0 s

mode **D**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	60 min.
<b>Hold</b>	3 s
<b>Interval</b>	3 s

mode **E**

- 
- 
- 
- 
- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	60 min.
<b>Hold</b>	1 s
<b>Interval</b>	2 s

mode **F**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	45 min.
<b>Hold</b>	2 s
<b>Interval</b>	2 s

mode **G**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	60 min.
<b>Hold</b>	2 s
<b>Interval</b>	2 s

mode **H**

- 
- 
- 
- 
- 
- 
- 
- 

<b>Time</b>	30 min.
<b>Hold</b>	2 s
<b>Interval</b>	3 s